

MDR[®] Video for Speech Pathology

The ideal recording solution for modified barium swallow studies

PACSGEAR, the leader in medical imaging connectivity has developed **MDR[®] Video for Speech Pathology**—designed with the imaging needs of Speech Pathologists in mind.

MDR Video for Speech Pathology provides the following features:

- ⚙ Synchronized audio and video can be recorded to CD/DVD/USB media and/or sent to PACS
- ⚙ Recording can be triggered from the fluoroscopic footswitch or an optional footswitch
- ⚙ Video clips can be trimmed and single frames captured so that only diagnostically relevant data is recorded to media or sent to PACS
- ⚙ **MDR Video RP** prints full color inkjet labels including patient and study data onto CDs/DVDs
- ⚙ Programmable function keys allow images to be easily labeled with text overlays to identify consistency of the swallowed substance
- ⚙ The **MDR Video Cart** allows for easy mobility in and out of the imaging room and review area
- ⚙ A touch screen monitor is provided for user interface and playback
- ⚙ Viewing software is recorded on every disk



 MDR[®] Video RP



 MDR[®] Video



MDR Video RP Cart

Get Started

For more information about **MDR Video for Speech Pathology**, contact us at +1 925 225 6100 or email getstarted@pacsgear.com