

 **PacsConnect**[™]

Pacs**Connect**[™] is a proven workflow solution for hospitals and imaging centers. Built on open standards such as DICOM and HL7, Pacs**Connect** enables IHE workflow for interoperability with your existing PACS. Simple and powerful software delivers data accurately and efficiently throughout your enterprise. Most importantly, PACS Administrators can monitor and configure the application anywhere on the network.

Pacs**Connect** lets technologists select the appropriate study at each modality without having to manually re-enter patient demographics, which eliminates data entry errors. Includes full DICOM and HL7 trace logs, worklist/patient queries, and lookup tables for easy data re-mapping.

Proven DICOM Modality Worklist (MWL) Technology

Why PacsConnect?

- ⚙ **Easy** – Simple user interface manages HL7 to DICOM MWL mappings
- ⚙ **Fast** – Quick, wizard-based installation and configuration
- ⚙ **“It just works”** – Validated against a wide variety of PACS and modalities




PacsConnect Features

- ⚙ Proven DICOM Modality Worklist (MWL) technology – validated against a wide variety of PACS and modalities
- ⚙ Quick, wizard-based installation
- ⚙ Rules-based HL7 forwarding with flexible field mappings
- ⚙ Supports all messages for HL7 versions 2.2, 2.3 and 2.4
- ⚙ Customized integration with non-HL7 RIS systems
- ⚙ Manage remotely from any PC

System Requirements

- ⚙ Microsoft® Windows 2000 or 2003 Server
- ⚙ Microsoft SQL Server
- ⚙ 1GB memory minimum
- ⚙ 30GB disk minimum (RAID recommended)
- ⚙ UPS



PacsConnect is easy to install and configure

PACSGEAR – the leader in Document and Multimedia Connectivity for PACS/EHR

Get Started

For more information about PacsConnect, contact us at +1 925 225 6100 or email getstarted@pacsgear.com