

 Print2PACS™

Print2PACS™ lets you convert and send documents and images from any Windows® 2000 or XP application that supports Windows-based printing. Eliminates paper, improves workflow and provides access to important clinical documentation. Now you can send Word files, PDFs, reports, email and EHR/EMR, cardiology and RIS information directly to your PACS.


Print2PACS is software that appears as a printer selection within your application. Select “Print2PACS” and print documents to PACS as you would print on paper. Cath labs, cardiology and DEXA/BMD scanners are just a few of the application areas where Print2PACS improves workflow. In PACS, view documents along with patient studies at the point of care.

As easy as clicking File > Print > Print2PACS

Why Print2PACS?

- ✦ **Easy** – No training required
- ✦ **Fast** – Send Windows-printable documents to PACS/EHR in less than one minute
- ✦ **“It just works”** – Converts a wide variety of documents and images. Send Word files, PDFs, reports, email, and EHR/EMR, cardiology and RIS information directly to your PACS.

How it works... As easy as 1, 2, 3!!!



1. From your Windows-based application, select **File > Print** and then select **"Print2PACS"** as the printer



2. Select patient/study demographics from your RIS or DICOM archive and click **"Continue"** to convert and send the document to PACS/EHR



3. The document is now matched with the appropriate study

Print2PACS Features

- ⚙ Supports color and grayscale imaging.
- ⚙ Interfaces to RIS worklist or DICOM Archive to provide automatic loading of patient demographics for ease of use. DICOM Modality Worklist (MWL) ensures an accurate match with patient studies.
- ⚙ Simple user interface allows anyone to add documents without prior knowledge of PACS. No training required.
- ⚙ Allows import of documents or images from any Windows 2000 or XP application that supports Windows-based printing.
- ⚙ Optional image invert feature reduces display glare.
- ⚙ Easily set the DPI to change the size and quality of all printed output sent to PACS/EHR.



Need electronic forms or document scanning? Consider **PacsSCAN™** – industry leading document scanning software.

Get Started

For more information about Print2PACS, contact us at +1 925 846 9600 or email getstarted@pacsgear.com