

 **Stacker**[™]

Stacker[™] automates the manual process of importing stacks of DICOM CDs/DVDs. Using reliable robotic disc readers, **Stacker** can import up to 50 discs in a single session. PACS Administrators and front office personnel alike will appreciate **Stacker**'s power and simplicity. Just load discs and click "Start."

While reading discs, **Stacker** readily identifies unreadable CDs/DVDs. **Stacker Pro** sends rejected CDs/DVDs to a separate "reject chute" for continuous, unattended operation. After reading, successfully read discs are displayed on a queue for final processing and can be quickly localized with patient information from DICOM Modality Worklist or DICOM archives and sent to PACS for comparison with current studies. For additional flexibility, send studies to **GEARView QC**, PACSGEAR's new Quality Control workstation.

Automates DICOM CD/DVD import

Why Stacker?

- ⚙ **Easy** – Imports up to 50 discs in a single session
- ⚙ **Fast** – Sends to PACS or GEARView QC for review / editing
- ⚙ **"It just works"** – Uses reliable Primera[®] Bravo[®] robotic systems

How Stacker Works



1. Load DICOM CDs and/or DVDs and click "Start."



2. Localize patient information using DICOM Modality Worklist and send to PACS or GEARView QC.

Stacker Features

- ⚙ Imports up to 50 discs in a single session.
- ⚙ Rejected discs are easily identified. Stacker Pro sends rejected discs to a separate location for continuous, unattended operation.
- ⚙ Successfully read discs are displayed on a queue for final processing, including media sent from other PacsSCAN locations.
- ⚙ Interface to RIS worklist or DICOM Archive to provide automatic loading of patient demographics for ease of use. DICOM Modality Worklist (MWL) ensures an accurate match with patient studies.
- ⚙ Includes PacsSCAN document scanning – scan and send documents to PACS/EHR using TWAIN compliant scanners.
- ⚙ Can be combined with GEARView QC for additional flexibility.

Features Reliable Primera® Bravo® Robotic Systems



⚙ Stacker™ Pro



⚙ Stacker™

Get Started

For more information about **Stacker**, contact us at +1 925 846 9600 or email getstarted@pacsgear.com